

***C. elegans* Slow-killing assay (SKA) Large Plate Selection Protocol**

by Erik Andersen June 22, 2010

Assay timing (each step is described below)

- Day 1: Streak out PA14 from frozen stock.
Day 2: Pour SKA plates, inoculate PA14 to LB (at 2 PM)
Day 3: Spot PA14 to center of SKA plate, spread to cover 75% of plate, and put plates at 37°C
Day 4: Take plates out of 37°C (at 2 PM)
Day 5: Put 1000-2000 L4s on each SKA plate
Days 9: Collect worms alive after 100 hours on SKA plates

Plate preparation

1. Prepare SKA plate medium and autoclave

	250 mL	500 mL	1 L
NaCl	0.75 g	1.5 g	3.0 g
BactoAgar	4.25 g	8.5 g	17 g
Peptone	0.875 g	1.75 g	3.5 g
Sterile water	243.75 mL	487.5 mL	975 mL

Autoclave, allow to cool to 55 °C and add the following

Cholesterol (5 mg/mL in EtOH)	0.25 mL	0.5 mL	1 mL
1 M CaCl ₂	0.25 mL	0.5 mL	1 mL
1 M MgSO ₄	0.25 mL	0.5 mL	1 mL
1 M KH ₂ PO ₄ (pH 6)	6.25 mL	12.5 mL	25 mL
FUDR, filter sterile (100 mg/mL)	125 µL	250 µL	500 µL

2. Pour 25 mL per 10 cm plate. Flame tops of plates to remove bubbles, if needed.
3. Inoculate 100 ml LB with *Pseudomonas ariginosa* (strain PA14) from a freshly streaked plate.
4. Grow at 37°C overnight for 24 hours.
5. Spot 800 µL of PA14 onto each plate, spread to cover 75% of plate, and put the plates at 37°C for 24 hours in a closed box. Be careful to not scratch the tops of the plates.
6. Remove the plates from 37°C and keep at room temp for 24 hours. Use plates immediately.

Slow-killing Assay Large Plate Selection

1. Put 1000-2000 L4 hermaphrodites onto each SKA plate.
2. Put the plates at 25°C.
3. At 96 hours, collect the living worms. How?
4. Isolate DNA by lysis
5. Make libraries
6. Run sequence
7. Look for enrichment of genomic regions in neighborhood of *npr-1*.