



Show document outline

# CRISPR-Cas9 genome editing by RNP injection

By Shanon Brady, Daehan Lee, Clay Dilks, Emily Koury, and Erik Andersen

Updated November, 2021

This protocol describes the steps to make a CRISPR-Cas9 mixture for microinjection and the steps after microinjection to screen for edits. Microinjection is a separate protocol. Please see the **ordering information guide** document for Andersen lab-specific ordering instructions for the reagents below. For a quick video explaining the basic process please follow [this link](#). Note that there are key differences in the video to this protocol and you should follow this protocol explicitly where they may differ.

## Reagents:

### sgRNA:

Design your guides near your targeted edit site using the instructions below. Order sgRNAs from Synthego and resuspend with provided TE to 100  $\mu$ M. Store at -20°C

### target repair oligo:

Design your repair oligo using the instructions below. Order as 4 nmol Ultramer DNA oligo from IDT and resuspend to 100  $\mu$ M in nuclease free water. Store at -20°C.

### *dpy-10* sgRNA:

*dpy-10* sgRNA sequence: GCUACCAUAGGCACCACGAG

Order from Synthego and resuspend with provided TE to 20  $\mu$ M. Aliquot and store at -20°C

### *dpy-10* repair oligo:

*dpy-10* repair oligo sequence:

CACTTGAACCTTCAATACGGCAAGATGAGAATGACTGGAAACCGTACCGCATGCGGTGCCTATGGT  
AGCGGAGCTTCACATGGCTTCAGACCAACAGCCTAT

Order using the above sequence as a 4 nmol Ultramer DNA oligo from IDT and resuspend to 10  $\mu$ M in nuclease free water. Aliquot and store at -20°C.

### Cas9:

40  $\mu$ M Cas9-NSL purified protein obtained from MacroLab Facility at the University of California, Berkeley. This order comes in 10  $\mu$ L aliquots that are marked with each use and discarded after four freeze/thaws.

